

Recommended Weight to Blood Draw Volume

[1lb:0.773mLs]

lbs	mLs
0.5	0.4
1	0.8
1.5	1.2
2	1.5
2.5	1.9
3	2.3
3.5	2.7
4	3.1
4.5	3.5
5	3.9
5.5	4.3
6	4.6
6.5	5.0
7	5.4
7.5	5.8
8	6.2
8.5	6.6
9	7.0
9.5	7.3
10	7.7
10.5	8.1
11	8.5
11.5	8.9
12	9.3
12.5	9.7
13	10.0
13.5	10.4
14	10.8
14.5	11.2
15	11.6
15.5	12.0
16	12.4
16.5	12.8
17	13.1
17.5	13.5
18	13.9
18.5	14.3
19	14.7
19.5	15.1

lbs	mLs
20	15.5
20.5	15.8
21	16.2
21.5	16.6
22	17.0
22.5	17.4
23	17.8
23.5	18.2
24	18.6
24.5	18.9
25	19.3
25.5	19.7
26	20.1
26.5	20.5
27	20.9
27.5	21.3
28	21.6
28.5	22.0
29	22.4
29.5	22.8
30	23.2
30.5	23.6
31	24.0
31.5	24.3
32	24.7
32.5	25.1
33	25.5
33.5	25.9
34	26.3
34.5	26.7
35	27.1
35.5	27.4
36	27.8
36.5	28.2
37	28.6
37.5	29.0
38	29.4
38.5	29.8
39	30.1

lbs	mLs
39.5	30.5
40	30.9
40.5	31.3
41	31.7
41.5	32.1
42	32.5
42.5	32.9
43	33.2
43.5	33.6
44	34.0
44.5	34.4
45	34.8
45.5	35.2
46	35.6
46.5	35.9
47	36.3
47.5	36.7
48	37.1
48.5	37.5
49	37.9
49.5	38.3
50	38.7
50.5	39.0
51	39.4
51.5	39.8
52	40.2
52.5	40.6
53	41.0
53.5	41.4
54	41.7
54.5	42.1
55	42.5
55.5	42.9
56	43.3
56.5	43.7
57	44.1
57.5	44.4
58	44.8
58.5	45.2

lbs	mLs
59	45.6
59.5	46.0
60	46.4
60.5	46.8
61	47.2
61.5	47.5
62	47.9
62.5	48.3
63	48.7
63.5	49.1
64	49.5
64.5	49.9
65	50.2
65.5	50.6
66	51.0
66.5	51.4
67	51.8
67.5	52.2
68	52.6
68.5	53.0
69	53.3
69.5	53.7
70	54.1
70.5	54.5
71	54.9
71.5	55.3
72	55.7
72.5	56.0
73	56.4
73.5	56.8
74	57.2
74.5	57.6
75	58.0
75.5	58.4
76	58.7
76.5	59.1
77	59.5
77.5	59.9
78	60.3